

Protein Analysis/Proteomics

ial and for use by approved personnel only. Copyright 2006 Cold Spring Harbor Laboratory
tribution. Do not copy without written permission from Cold Spring Harbor Laboratory Press.

The following section covers buffers used in protein analysis and investigation of protein–protein interactions. This section is presented in two parts: (1) basic methods in protein analysis for proteomics and (2) advanced protocols in identifying protein–protein interactions for proteomics. Areas covered include peptide sequencing, analysis of peptide/DNA complexes, and two-hybrid and library screening.

Caution: See Cautions Appendix for appropriate handling of materials marked with ▼.

Basic Methods in Protein Analysis for Proteomics

Chromatography

The following sections include solutions used in various chromatographic purification protocols.

■ *Size Exclusion Chromatography*

Denaturing Buffer Containing Urea PUR

8 M urea ▼

5 mM dithiothreitol ▼

150 mM NaCl

50 mM Tris-Cl (pH 7.5)

Warm the solution to 40°C to dissolve the urea. Adjust the pH to 7.5 if necessary.

Denaturing Buffer Containing Guanidine Hydrochloride PUR

6 M guanidine hydrochloride ▼
 5 mM dithiothreitol ▼
 50 mM Tris-Cl (pH 7.5)

Warm the solution to 40°C to dissolve the guanidine hydrochloride. Adjust the pH to 7.5 if necessary.

Material and for use by approved personnel only. Copyright 2006 Cold Spring Harbor Laboratory Press. Do not copy without written permission from Cold Spring Harbor Laboratory Press.

Peptide Elution Buffers PUR

20% acetonitrile ▼ / 0.1% TFA ▼ in H₂O
 0.1% TFA in H₂O
 or
 0.1% acetic acid ▼ in H₂O

■ **Hydrophobic Interaction Chromatography****High-density Phenyl-substituted Matrix Column** PUR

Suitable columns include the POROS HP2 (M) chromatography column (Applied Biosystems), SOURCE 15 PHE PE (100 × 4.6 mm), or Bio-Gel TSK Phenyl-5PW HIC column (75 × 7.5 mm I.D. 1000 Å) (Bio-Rad).

Dimethyl Pimelimidate (DMP) Solution PUR

Prepare 10 ml of 20 mM DMP (Pierce) in 0.1 M triethanolamine ▼ (pH 9.5; the pH of the solution decreases to ~8.0 due to the acidity of DMP). This reagent must be prepared fresh immediately before use.

Column Regeneration Buffer PUR

10 mM Tris-HCl (pH 7.6)
 1 mM EDTA
 2.5 M NaCl
 1% Nonidet P-40 (NP-40)

NP-40 will never completely dissolve. Swirl just before use.

Column Storage Buffer PUR

10 mM Tris-HCl (pH 7.6)

1 mM EDTA

0.3 M NaCl

0.04% sodium azide ▼

ial and for use by approved personnel only. Copyright 2006, Cold Spring Harbor Laboratory.
 distribution. Do not copy without written permission from Cold Spring Harbor Laboratory Press.

HEMGN Buffer PUR25 mM HEPES (K⁺, pH 7.6)

0.1 mM EDTA

2.5 mM MgCl₂ ▼

10% glycerol

0.025% Nonidet P-40

1 mM dithiothreitol (DTT) ▼

Add DTT fresh immediately before use. Store the buffer at 4°C.

■ **Immobilized Metal-ion Affinity Chromatography****Supplemented Buffer P** PUR

50 mM sodium phosphate (pH 7.4) ▼

0.5 M NaCl

0.5 mg/ml lysozyme ▼

1 mM PMSF ▼

1.7 units/ml Benzonase (Merck)

1 mM MgCl₂ ▼

Add the lysozyme, PMSF, and Benzonase just before use.

Binding Buffer PUR

20 mM sodium phosphate (pH 7.4) ▼

0.5 M NaCl

20 mM imidazole ▼

Elution Buffer PUR

20 mM sodium phosphate (pH 7.4) ▼

0.5 M NaCl

500 mM imidazole ▼

Inclusion Body Wash Buffer PUR

2 M urea ▼
 500 mM NaCl
 2% Triton X-100 ▼
 20 mM Tris-HCl (pH 8.0)

ial and for use by approved personnel only. Copyright 2006 Cold Spring Harbor Laboratory
 istribution. Do not distribute without written permission from Cold Spring Harbor Laboratory Press.

20 mM imidazole ▼
 500 mM NaCl
 1 mM 2-mercaptoethanol ▼
 20 mM Tris-HCl (pH 8.0)

Urea Binding Buffer PUR

For preparation of buffers containing urea, weigh the denaturant, add 2 M imidazole ▼, 1 M Tris-HCl (pH 8.0), and H₂O to 90% of the final volume. Heat the container, which is cooled by the solubilization process of the denaturant, to room temperature in a water bath at 45°C while stirring. Avoid high temperatures. The denaturant will dissolve within minutes. Alternatively, stir overnight at room temperature. Adjust the pH, add H₂O to the final volume, and filter the solution. The pH adjustment must be done at room temperature. Add 2-mercaptoethanol ▼ immediately before use.

■ **Chromatofocusing**

The following section includes a solution used in purifying proteins by chromatofocusing.

Acid Solutions PUR

2 M acetic acid ▼
 2 M HCl ▼
 iminodiacetic acid (saturated solution) ▼

These solutions are used for adjusting the pH of start and elution buffers.

■ Two-dimensional Gel Electrophoresis

The following section includes solutions for preparation and visualization of proteins by two-dimensional (2D) electrophoresis.

Thiourea/Urea Lysis Buffer PUR

ial and for use by approved personnel only. Copyright 2006 Cold Spring Harbor Laboratory
 ution. Do not copy without written permission from Cold Spring Harbor Laboratory Press.

- 2 M thiourea ▼
- 7 M urea ▼
- 4% (w/v) CHAPS ▼
- 1% (w/v) DTT ▼
- 2% (v/v) carrier ampholytes (pH 3–10)
- 10 mM Pefabloc proteinase inhibitor

To prepare 50 ml of thiourea/urea lysis buffer, dissolve 22.0 g of urea (GE Healthcare) in deionized H₂O, add 8.0 g of thiourea (Sigma-Aldrich), and adjust the volume to 50 ml with deionized H₂O. Add 0.5 g of Serdolit MB-1 mixed-bed ion-exchange resin (Serva), stir for 10 minutes, and filter. To 48 ml of the urea solution add 2.0 g of CHAPS, 1.0 ml of Pharmalyte (pH 3–10; GE Healthcare), 0.5 g of DTT, and, immediately before use, 50 mg of Pefabloc proteinase inhibitor (VWR).

Urea Lysis Buffer PUR

- 9.5 M urea ▼
- 1% (w/v) DTT ▼
- 2% (w/v) CHAPS ▼
- 2% (v/v) carrier ampholytes (pH 3–10)
- 10 mM Pefabloc proteinase inhibitor

To prepare 50 ml of urea lysis buffer, dissolve 30.0 g of urea (GE Healthcare) in deionized H₂O and adjust the volume to 50 ml. Add 0.5 g of Serdolit MB-1, stir for 10 minutes, and filter. Add 1.0 g of CHAPS, 0.5 g of DTT, and 1.0 ml of Pharmalyte (pH 3–10) to 48 ml of the filtered urea solution. If necessary, add 50 mg of Pefabloc proteinase inhibitor immediately before use.

IMPORTANT: Lysis buffer must be freshly prepared. Alternatively, make 1-ml aliquots and store at –80°C for up to several months. Do not refreeze lysis buffer that has already been thawed once! Never heat urea solutions above 37°C. Otherwise, protein carbamylation may occur.